

Silencing cargo gene expression in lentiviral vector and AAV production to improve vector yield and purity

Qian Liu¹, Sherin Parokkaran Johny¹, Matthew Tridgett, Mohine Alam, Fran Bennett, Siobhan Clerkin, Michal Fleszar, Meenakshi Raghunath¹, Maria Patricio, Michael White, Weiheng Su², Vaughan Leydon¹


¹Minaris Advanced Therapies, London, UK, ²Independent Consultant, Oxford, UK.

Poster


1. Introduction

Lentiviral vectors (LVVs) and adeno-associated viral vectors (AAVs) are widely used as delivery platforms for genetic modification in both *ex vivo* cell therapies and *in vivo* gene therapies. However, their manufacturing can be complicated by the expression of certain therapeutic genes, particularly when these genes are cytotoxic. In such cases, cell death during production can generate byproducts that increase the burden on downstream processing. These challenges are difficult to avoid because cargo gene expression is typically high, driven by the large quantities of transfer plasmid used during manufacturing. In addition, strong viral promoters positioned upstream of the cargo gene can unintentionally amplify expression within the production system. As a result, there is a need to develop strategies to suppress cargo gene expression during production, enabling reliable manufacturing of viral vectors that carry potent but cytotoxic therapeutic genes.


2. Mechanism of LVV cargo gene silencing



3. LVV Cargo gene knockdown improves production cell health & reduces DSP burden




4. LVV Cargo gene silencing increases titres




5. LVV Cargo gene knockdown facilitates stable producer cell line development

- Constitutive expression of certain cargo genes may limit producer cell growth so a stable cell line cannot be developed.
- We stably integrated either conventional or silenced anti-CD19 CAR transfer plasmids into XOFLX[®] Packaging cells.
- Silencing the cargo gene enabled generation of stable producer cells that otherwise were unrecoverable (Figure 4 A-C).




6. AAV Gene Silencing Technology against EGFP reporter

- shRNA approaches to gene silencing in AAV may lead to PKR activation and translational shutdown.
- Our novel AAV gene silencing technology works with the biology of AAV leading to demonstrable knockdown of reporter genes by fluorescence microscopy (Figure 5A) and flow cytometry (Figure 5B) without affecting AAV yields (Figure 5C).



7. AAV Gene Silencing of Bax during plasmid transfection

- Different variations of the AAV silencing strategy were assessed using a known toxic transgene, BAX (BCL2 associated X, apoptosis regulator).
- Modest improvements (~1.4x) in cell viability during AAV production process (Figure 6B) were observed from the two indirect targeting approaches.
- A 16-fold improvement in AAV yield was observed for one strategy (Figure 6C).
- This technology has the potential to enable AAV manufacture for previously intractable human transgenes where shRNA-based approaches have been unsuccessful.



8. Conclusions & Future work

- Cargo gene expression levels were reduced significantly during lentiviral and AAV vector production processes, using different silencing mechanisms developed by Minaris.
- LVV and AAV production yields were maintained (for non-toxic genes) or increased (for cytotoxic transgenes) by implementing the corresponding cargo gene silencing technology.
- Production cell viability was improved, reducing cell death-related process impurities and potentially enabling more efficient and robust downstream process performance.
- The technologies can also facilitate the development of stable producer cell lines for viral vectors encoding cytotoxic cargo genes.
- We will test the silencing mechanisms on a wider range of transgenes in both transient and stable productions.
- Use of these cargo gene silencing systems will de-risk development of novel LVV- and AAV-based therapeutics, by improving production yield and reducing the unwanted protein products from the production process, resulting in decreased manufacturing costs and complexity and improved patient access to LVV- and AAV-based therapies.
- Cargo silencing systems are available for evaluation. Request at www.minaris.com

Corresponding author:
Qian.Liu@minaris.com

For information about our technologies,
please contact us at www.minaris.com

Scan the QR Code and
download the poster

OR CONTACT US

minaris.com/contact

